

2

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	

Snap Neck ND11

- ♦ The vials are used on instruments of Agilent, CTC, Thermo Scientific, Waters, etc.
- ♦ Snap Ring finish eliminates the need for crimping or decapping.
- ♦ The intergrated Micro-insert are available for in clear and amber glass, micro-inserts can be delivered pre-assembled in vials.
- ♦ Snap ring caps are also available in a soft quality with transparent or blue caps, soft caps are more convenient in handling (pushing on/removal).
- ♦ The septa shapes have Single and Snow type to choose.



Ref.No.	C0000020	C0000021	C0000022	C0000023	C0000130
Description	2mL Clear Glass Flat Base 11mm Snap Vial Wide Opening.	2mL Clear Glass Flat Base 11mm Snap Vial Wide Opening with Label.	2mL Amber Glass Flat Base 11mm Snap Vial Wide Opening.	2mL Amber Glass Flat Base 11mm Snap Vial Wide Opening with Label.	2mL 11mm Snap Top Clear Vial Bottom with Integrated 0.2ml Glass Micro-insert.
100pcs/pk					



Ref.No.	C0001023 ★	C0001827 ★	C0001867/C0001868 ★	C0001870 ★
Description	2mL 11mm PP Snap Top Vial with 0.3mL Micro-Vial Transparent. 12x32mm.	2mL 11mm PP Snap Top Vial with 0.3mL Micro-Vial Amber. 12x32mm.	450uL Filter Vial with 0.45um/0.22um PTFE(Hydrophilic) membrane and a CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	2mL 11mm PP Snap Top Vial with 0.7mL Micro-Vial Transparent. 12x32mm.
100pcs/pk				

>ND 11 Closures



Ref.No.	C0000170	C0000171	C0000172	C0000173
Description	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	BLUE 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	RED 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone/Red PTFE Septa 1mm Thick.
100pcs/pk				



Ref.No.	C0000174	C0000175	C0000176	C0000177
Description	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick PRE-SLIT.	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick PRE-SLIT (snow shape).	CLEAR 11mm Open Top Snap Cap with Transparent TEF/Natural Rubber red-orange Septa 1mm Thick.	CLEAR 11mm Open Top Snap Cap with White PTFE/Red Silicone Septa 1mm Thick.
100pcs/pk				

11mm Snap Cap With Septa

PTFE/Silicone septa are most popular for HPLC applications.



Product Details

Zhejiang ALWSCI Technologies Co.,Ltd is one of the leading China 11mm snap cap with septa manufacturers and suppliers, welcome to wholesale cheap ptfе seal, teflon septa from us.

FEATURES

- 1.PTFE/Silicone septa are most popular for HPLC applications, and Pre-slit septa are easier to pierce with needles.
- 2.Usually, the cap color is clear or blue, other special colors can be also produced accordingly.
- 3.Easy to apply and easy to remove.
- 4.Pre-assembled caps and septa are convenient.

P/N	Size	Packing	Description
11-SP4001	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with Red PTFE/White Silicone 1mm (.040") thick
11-SP4001-B	11mm Snap	100/PK	Blue,11mm Open Top Snap Cap with Red PTFE/White Silicone 1mm (.040") thick
11-SP4001-R	11mm Snap	100/PK	RED, 11mm Open Top Snap Cap with Red PTFE/White Silicone 1mm (.040") thick
11-SP4002	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with Red PTFE/White Silicone/Red PTFE 1mm (.040") thick
11-SP4003	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with Red PTFE/White Silicone 1mm (.040") thick, PRE-SLIT
11-SP4004	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with Red PTFE/White Silicone 1mm (.040") thick, PRE-SLIT, Snow Shape
11-SP4005	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with Transperant PTFE/Natural Rubber red-orange Septa 1mm (.040") thick
11-SP4006	11mm Snap	100/PK	CLEAR, 11mm Open Top Snap Cap with White PTFE/Red Silicone 1mm (.040") thick

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	