

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	

Crimp Top Headspace Vials gas chromatography mass spectrometry

Crimp headspace vial widely used in Headspace and gas chromatography mass spectrometry analysis.

Product Details

1) Description

10ml and 20ml crimp headspace vial is widely used in Headspace and GC, GCMS application.

Item No.	Description
C0000033	10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial. 100pcs/pk.
C0000034	10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Rounded-flat Bottom. Crimp Headspace Vial. 100pcs/pk.
C0000035	10mL Amber Glass 22.5x46mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial. 100pcs/pk.
C0000036	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial. 100pcs/pk.
C0000038	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. 100pcs/pk.
C0000039	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. Short Neck. 100pcs/pk.
C0000040	20mL Amber Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. Short Neck. 100pcs/pk.
C0001877	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial with Graduations. 100pcs/pk.

2) Details

Product Name Headspace vial

Volume 10 ml and 20 ml

Neck type Crimp(ND20mm)

Vial Color Clear& Amber

Cap Material Aluminum or Bi-metallic

MOQ 1000pcs

Packing 100pcs/pack

Delivery 1-2 weeks

Certificate ISO9001,SGS,COA

20mm Crimp Pressure Release Cap

20mm Open Top Aluminum Crimp Cap (pressure release) with Septa



Product Details

1. Product Description

- 20.2mm Open Top Aluminum Crimp Cap (pressure release) with Natural PTFE/White Septa. Used for normal volatile substance, with excellent reseal ability and wide chemical compatibility.
- Pressure release seals (Aluminum-PR) are designed to open when the internal pressure reaches a certain value.
- Different Septa for customer choice. Suitable for Crimp Headspace vials for the analysis of volatile organic solids and liquid samples after gasification.

2. Specification of Crimp caps with PTFE/Silicone Septa

P/N	Packing	Description
C0000315	100pcs/pk	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick.
C0001296	100pcs/pk	20.2mm Open Top Aluminum Crimp Cap (pressure release) with Natural PTFE/White Septa 3mm Thick. 100pcs/pk.
C0001389	100pcs/pk	20.2mm Open Top Aluminum Crimp Cap (pressure release) with 20mm Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick.
C0001403	100pcs/pk	20.2mm Open Top Aluminum Crimp Cap (pressure release) .

Teflon Faced Septa

ND20 Crimp Cap with teflon faced septa is used for Crimp headspace vial. It can be used in autosampler like CTC, PerkinElme.



Product Details

1) Description

ND20 Crimp Cap with teflon faced septa is used for Crimp headspace vial. It can be used in autosampler like CTC, PerkinElme.

Ref.No.	Description
C0000196	20mm Aluminum Cap with Natural PTFE/White Silicone Septa
C0000198	20mm Aluminum Cap with White PTFE/White Silicone Septa
C0000205	20mm Gold Magnetic Crimp Cap with Natural PTFE/White Silicone Septa
C0000206	20mm Bi-metallic Crimp Cap Blue & Silver with Natural PTFE/White Silicone Septa
C0000835	20mm Gold Magnetic Crimp Cap with White PTFE/Blue Silicone Septa
C0000951	20mm Bi-metallic Crimp Cap Red & Silver with Natural PTFE/Natural Silicone Septa
C0000952	20mm Bi-metallic Crimp Cap Gold & Silver with Natural PTFE/Natural Silicone Septa

2) Details

Cap Material	Aluminum/Metal
Cap Size	20mm
Septa Material	PTFE/Silicone
Application	Used for ND20 Crimp Headspace Vial
Packing	100pcs/pack

3) Advantages

- * High chemical resistance
- * Butyl septa are available
- * Colors can be customized
- * Other diameter are available

Headspace Aluminum Crimp Caps

20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.



Product Details

1. Product Description

- Aluminum Crimp Cap with PTFE/Silicone septa. Used for normal volatile substance, with excellent reseal ability and wide chemical compatibility.
- Different colors for customer choice. Suitable for Crimp Headspace vials for the analysis of volatile organic solids and liquid samples after gasification.



2. Specification of Crimp caps with PTFE/Silicone Septa

P/N	Packing	Description
C0000196	100pcs/pk	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick.
C0000197	100pcs/pk	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Silicone Septa 3mm Thick.
C0000198	100pcs/pk	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/White Silicone Septa 3mm Thick.
C0000199	100pcs/pk	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/Blue Silicone Septa 3mm Thick.

Aluminum Crimp Caps For GC



20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.

Product Details

1. Product Code & Description

20mm Crimp seals with prefitted septa for headspace vials Seals.

Item Number	Description
C0000198	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.

* Product Introduction

Constructed of aluminum, these crimp vials are ideal for use with your crimp neck vials. Aluminum crimps tightly around the neck for a tight seal.

2. Details

Product Name	Aluminum Crimp Caps
Size	20mm
Material	Aluminum, Magnetic, Bi-Metallic , PTFE/Silicone
Application	GC, Crimp Headspace vial
Feature	Autoclavable, easy sterilization
MOQ	1000pcs
Packing	100pcs/pack, Inner PE bag; Standard Carton outside

3. Advantages

1. Caps are ready to use with any headspace crimp top vials.
2. Pre-assembled caps and septa are convenient and minimize contamination from handling.
3. It is autoclavable, easy sterilization.
4. Very low background peaks, very low permeability.
5. Both Aluminum, Magnetic, Bi-Metallic caps are available.

20mm Crimp Caps GC MS

20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.



Product Details

1. Production details

GC Headspace technology for the analysis of volatile organic solids and liquid samples after gasification.

Item Description:

Item Number	Description
C0000196	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.
Product Name	20mm Crimp Cap
Size	Diameter: 20mm
Material	Aluminum, PTFE/Silicone
Application	Used for headspace crimp vial for GC analysis in laboratory
Feature	Various color 's septa for option
MOQ	1000pcs
Packing	100pcs/bag, Inner PE bag; Outer: Standard Carton

Advantage

1. Crimp headspace Aluminum cap, suitable for GCMS.
2. Types of liner, like Butyl/PTFE,Pharma-fix liner, PTFE/Silicone, Butyl stopper.
3. Produced by ALWSCI, it have high quality control.
4. Provide OEM,ODM services.

20mm Ultra Clean Series

Ultra Clean 20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Silicone Septa 3mm Thick.



Product Details

1. Product Details

UltraClean PTFE/silicone septa is optimized for gas chromatography and mass spectrometry.

Item number	Description
C0000973	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Transparent Silicone Septa 3mm Thick (UltraClean). 100pcs/pk.



2. Details

Product name	20mm Ultra Clean Crimp Cap
Size	20mm
Material	Aluminum Cap, PTFE/Silicone Septa
Application	gas chromatography and mass spectrometry
Features	UltraClean. Various Color Choices
MOQ	1000pcs
Packing	100pcs/pack Caps and septa packed in transparent PE bag. Neutral carton outside, pallet is available to protect quality better. OEM packing can be also provided.

Crimp Top Headspace Cap & Septa

20.2mm Open Top Aluminum Crimp Cap (10mm hole) with Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick. 100pcs/pk



Product Details

1. Product Description

- Aluminum Crimp Cap with Butyl Septa. Effectively withstand most solvents and reduce dross, can withstand multiple penetration.
- The gray butyl outer sealing edge is tightly seamed with the sample bottle mouth to prevent the loss of fixed gas and the sealing is better.
- Suitable for Crimp Headspace vials for the analysis of volatile organic solids and liquid samples after gasification.

2. Specification of Crimp caps with butyl septa

P/N	Packing	Description
C0000855	100pcs/pk	20.2mm Open Top Aluminum Crimp Cap (10mm hole) with Moulded Grey PTFE/Butyl Septa one ring 3mm Thick.
C0001012	100pcs/pk	20.2mm Open Top Aluminum Crimp Cap (10mm hole) with Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick.

3. Specification of Butyl Septa

P/N	Packing	Description
C0000407	100pcs/pk	20mm Moulded Grey PTFE/Butyl Septa one ring 3mm thick.
C0000448	100pcs/pk	20mm Moulded Grey Butyl Septa one ring. 3mm thick.100pcs/pk.
C0000449	100pcs/pk	20mm Moulded Grey Butyl Septa both rings. 3mm thick.
C0000450	100pcs/pk	20mm Grey PTFE/Butyl without ring. 3mm thick.
C0000423	100pcs/pk	20mm Moulded Pharma-fix-septa Grey PTFE/Butyl 3mm thick.



Open Top Gold Magnetic Crimp Cap

20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.

Product Details

ALWSCI: Leading Manufacturer in HPLC/GC Consumables

We are a factory, not trading company

Quick Details

Classification: Laboratory Bottl

Model Number: C0000205

P/N: 20-SP3004-GM

Diameter: 20mm

Hole: 10mm

Vial Top: Crimp

Sample: Provided as free

Ref.No.	P/N	Description
C0000205	20-SP3004-GM	20mm Open Top Gold Magnetic Crimp Cap (10r th 20mm Natural PTFE/White Silicone Septa 3mm cs/pk.

20mm Open Top Crimp Cap GC MS

20mm open top crimp cap is perfectly match with 6ml/10ml/20ml headspace vial.



Product Details

Product Description:

Item No.	Description
C0000196	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.
C0000199	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm White PTFE/Blue Silicone Septa 3mm Thick. 100pcs/pk.
C0000206	20mm Open Top Bi-Metallic Crimp Cap (8mm hole) Blue & Silver with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.
C0000205	20mm Open Top Gold Magnetic Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.

Product Specification:

Product Name	20mm Open Top Crimp Cap GC MS
Cap choices	Aluminum or metal
Septa choices	Ptfe/silicone
Thread	Crimp ND20mm
Used	Used to 6ml/10ml/20ml headspace vial
Package	100PCS/PE Bag; 100Bags/CTN

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Казань (843)206-01-48	Новокузнецк (3843)20-46-81	Смоленск (4812)29-41-54
Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Новосибирск (383)227-86-73	Сочи (862)225-72-31
Астрахань (8512)99-46-04	Калуга (4842)92-23-67	Омск (3812)21-46-40	Ставрополь (8652)20-65-13
Барнаул (3852)73-04-60	Кемерово (3842)65-04-62	Орел (4862)44-53-42	Сургут (3462)77-98-35
Белгород (4722)40-23-64	Киров (8332)68-02-04	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Брянск (4832)59-03-52	Краснодар (861)203-40-90	Пенза (8412)22-31-16	Томск (3822)98-41-53
Владивосток (423)249-28-31	Красноярск (391)204-63-61	Пермь (342)205-81-47	Тула (4872)74-02-29
Волгоград (844)278-03-48	Курск (4712)77-13-04	Ростов-на-Дону (863)308-18-15	Тюмень (3452)66-21-18
Вологда (8172)26-41-59	Липецк (4742)52-20-81	Рязань (4912)46-61-64	Ульяновск (8422)24-23-59
Воронеж (473)204-51-73	Магнитогорск (3519)55-03-13	Самара (846)206-03-16	Уфа (347)229-48-12
Екатеринбург (343)384-55-89	Москва (495)268-04-70	Санкт-Петербург (812)309-46-40	Хабаровск (4212)92-98-04
Иваново (4932)77-34-06	Мурманск (8152)59-64-93	Саратов (845)249-38-78	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Набережные Челны (8552)20-53-41	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Иркутск (395)279-98-46	Нижний Новгород (831)429-08-12	Симферополь (3652)67-13-56	Ярославль (4852)69-52-93
Россия (495)268-04-70	Киргизия (996)312-96-26-47	Казахстан (7172)727-132	